

Whole Mount Mouse Embryo in situ Protocol

Day 0

1. Dissect embryos in PBS
2. Fix embryos in 4% paraformaldehyde at room temperature or 4C.
3. Dehydrate the embryos into 100% methanol and store at -20C. (optional)
4. (optional) Prior to use, bleach the embryos in 5:1 methanol: 30% hydrogen peroxide or 5:1 PFA:30% hydrogen peroxide for 5-6 hours at room temperature, followed by several washes in methanol or PFA. The embryos can then be stored in methanol at -20C, or in PFA at 4C indefinitely.

Day 1

- Each wash, unless indicated, is with at least 2ml of solution
 - Use an inverter mixer for all washes
1. If dehydrated, rehydrate the embryos through graded methanol series
 - a. 75% x 5min @4C
 - b. 50% x 5min @4C
 - c. 25% x 15min @4C (or until embryos sink*)
 - d. transfer samples to a 2ml eppendorf tube (Rnase, Dnase free microfuge tube)
 2. Wash three times in PBS
 - a. PBS x 5min
 - b. PBS x 5min
 - c. PBS x 5min
 3. Digest with 10ug/ml (2µl of 10mg/ml stock) of Proteinase K x 30 min @ RT*
 - a. Actual PK digestion time depends upon the age of the embryo. The following is a rough estimate of times:

E8.5 = 6 min
E9.5 = 10 min
E10.5 = 13 min
E11.5 = 15 min
E12.5 = 17 min
E14.5 = 20 min
E16.5 = 30 min
E18.5 = 40 min

*Deproteination is a critical step. Under-digestion will result in poor probe penetration. Over-digestion will result in a loss of structural integrity.

4. Wash in PBS
 - a. PBS x 1min
 - b. PBS x 5min
 - c. PBS x 5min
 - d. PBS x 5min
5. Discard wash solution*

*Pay special attention to discard as much of the PBT solution as possible.
6. Add 1.8 ml of pre-hybridization mix and allow the samples to incubate at the correct hybridization temp in a rotating hybridization oven for at least 1 hour.
 - a. As one hour nears, prepare the salmon sperm DNA and riboprobe as follows:
 - prepare 200ul of salmon sperm DNA per 2ml sample; denature at 90C for 10 min and place on ice until ready to use.
7. After at least 1 hour of prehybridization, add 200ul of denatured ssDNA to each sample, and about 100ng riboprobe to each sample. Hybridize overnight.

Day 2

1. Remove hybridization mix, replace with 2.0 ml 2X SSC x 10min @ hyb temp.
 - a. Wash with 2X SSC x 10min @ hyb temp
 - b. Wash with 2X SSC x 10min @ hyb temp
 - c. Wash with 2X SSC x 10min @ hyb temp
 - d. Incubate in 2X SSC x 60min @ 70C*

*70C kills any endogenous alkaline phosphatase activity. This is perhaps the most significant step towards reducing background. Breaking the 70C wash into two 30-minute sections may minimize damage to the tissue.
2. Wash with RNase A Buffer x 5 min @ RT
3. Discard, add RNase A Buffer + 2.0 ul of RNase A Enzyme (5mg/ml, 83U/mg) x 60min @37C
4. Discard, Wash four times with 1X Wash (BM)
 - a. 1X wash solution, x 10min @RT
 - b. 1X wash solution, x 10min @RT
 - c. 1X wash solution, x 60min @70C
5. Incubate in block buffer on inverter x 1 hr @ RT

- a. At this time, prepare 2.0ml of block buffer per sample and add 1 μ l Anti-digoxigenin antibody per sample (1:2000) to this master mix, invert briefly, and let sit at 4C until use.
6. After 1 hour, discard block buffer and add 2.0ml of pre-absorbed block buffer + antibody to samples and incubate overnight @ RT on an inverter.

DAY 3

1. Discard block solution.
2. Wash with 1X wash buffer and wash every hour throughout the day, for a minimum of 5 hours.
3. Wash with 1X was buffer overnight at RT.

DAY 4

1. Rinse with detection buffer for 10 min @ RT.
2. Detect with BM Purple at RT until desired signal strength is obtained (with longer probes you may be able to detect overnight).
3. Store Samples in 4% PFA at 4C.

Solutions

RNase Buffer:

0.5M NaCl (29mg/ml)

10mM Tris-HCl (pH 7.5)

Hybridization Mix:

50% Formamide by volume

50% 2 X SSC by volume

6% dextran sulfate by mass

SSC = saline – sodium citrate